Precautions

1. Clean exposed parts with a soft, slightly damp cloth. To remove stains, use a mild soap. Never use detergents, excess water, treated cloths, harsh cleaning agents, or sprays. Do not attempt to lubricate or open scale casing as this will void your warranty.
2. Although your scale is designed to be maintenance-free, it contains sensitive electronic parts. Avoid rough treatment. Do not drop, kick, or jump on it. Treat it with care to ensure the best performance.
3. This scale is designed to read weights between 6.6 and 330 lbs / 3 and 150 kgs / 6 ½ and 23 stone 8 lbs. To prevent damage, do not place more than 330 lbs / 150 kgs / 23 stone 8 lbs on it.
4. Do not store scale where you store cleaning chemicals. The vapors from some household products may affect the electronic components of your scale. Do not store scale on its side.
5. This scale is an extremely sensitive weighing device. In order to prevent run down of the battery, do not store anything on the scale.
6. Do not dispose of batteries in fire. Battery may explode or leak. Remove the batteries if the scale will not be used for a long period of time.

Caution: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Five (5) Year Limited Warranty

This scale is warranted against defects in materials of workmanship (excluding batteries) for five (5) years for the original purchaser from date of retail purchase. It does not cover damages or wear resulting from accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair. Do not return to retailer. For product support & warranty information: Please contact us with proof of purchase (purchase receipt) at:

Taylor Precision Products, Inc.
2220 Entrada Del Sol, Suite A
Las Cruces, New Mexico 88001 USA
Phone: 1-866-843-3905
Email: prodsupport@taylorusa.com

There are no express warranties except as listed above. This warranty gives you specific legal rights, and you may have other rights which vary from state to state. For additional product information, or warranty information outside the USA, please contact us through www.taylorusa.com. ©2015 Taylor Precision Products, Inc. and its affiliated companies. All rights reserved. Taylor® and Leading the Way in Accuracy® are registered trademarks of Taylor Precision Products, Inc. and its affiliated companies. All rights reserved.

Made to our exact specifications in China.
Not legal for trade.
www.taylorusa.com

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Electronic Scale

Instruction Manual

Thank you for purchasing a Taylor® scale. We strive to provide you with a high quality bath scale that meets your needs. Your opinion and experience with our product is important to us. If you are pleased with your Taylor product, please take a moment to share your thoughts with others by writing an on-line review. If you have questions or are not satisfied with your purchase, please contact us at prodsupport@taylorusa.com or call 1-866-843-3905.

Please read this instruction manual carefully before use. Keep these instructions handy for future reference.

General Operation

1. There is a strip under the battery cover to prevent battery drain in shipping. Please remove the strip before continuing with these instructions. Also, some models have a static cling label on the lens to prevent scratching. Please remove before use.
2. Your scale has been set at the factory to read weight in pounds. You may select kilograms or stone (1 stone = 14 pounds) units of measurement by pressing the button on the bottom of the scale.
3. Place scale on flat, hard surface. Carpeted or uneven floors may affect accuracy.
4. Step on the scale (disregarding the reading) and step off. This process insures the scale is set and is now ready for use. Do this whenever the scale is moved. **At all other times step straight onto the scale.**
5. Step onto the scale platform and remain still. The scale will count up to your weight. The display will flash twice when the weight is locked. Your weight is now displayed.
6. The scale will turn off automatically. NOTE: If you've weighed yourself and a "C" is displayed after the zero, weigh yourself again for a more accurate reading. Your weight cannot be locked-in if you do not stand still.

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
Battery Replacement

- Use 1- 3 volt CR2032 lithium battery (as installed).
- Remove screw (or screws) from the battery compartment cover on the bottom of the scale and remove the cover. Take care when removing and replacing the battery door screw. Do not turn the screw in the wrong direction or over tighten it as this may strip the screw threads in the plastic battery door.
- Remove the old battery from the compartment.
- Place the new battery with the "+" sign up into the battery compartment.
- Replace the battery compartment cover and screw(s).

NOTE: Please recycle or dispose of batteries per local regulations.

WARNING: Batteries may pose a choking hazard. As with all small items, do not let children handle batteries. If swallowed, seek medical attention immediately.

PRECAUTION: Do not dispose of batteries in fire. Batteries may explode or leak. Remove the battery if the scale will not be used for a long period of time.

Low Battery

Your Taylor® scale is equipped with a low battery indicator. It will display "Lo" or readings will grow dim or become irregular if the battery needs replacement.

Trouble Shooting

In Case of Difficulty...

- Place your feet in the same position every time you weigh yourself. Stand still during measurement.
- Weigh yourself without clothing whenever possible. Clothes can vary in weight.
- Move the scale to a location where it won't be bumped. If your scale is bumped while being activated or while in use, an inaccurate reading may occur.
- For the most consistent readings, weigh yourself at the same time of day. Allow for weight fluctuations due to food or drink in your system.
- "O_Ld" = overload warning.** The maximum capacity of this scale is 330 lbs x 0.2 lb / 150 kgs x 0.1 kg / 23 st 8 lb x ¼ lb. Should "O_Ld" appear on the display during a weighing, this indicates that the maximum weighing capacity of the scale has been exceeded. Remove the weight immediately; otherwise, permanent damage may occur.
- "Err" = Instability Error.** Ensure the scale is on a flat, stable surface. Step off and back onto the scale, standing still while your weight computes.

Reaching Target Weight

Maintaining proper weight is essential to healthy living. Obesity has been linked to a higher risk of heart disease, diabetes and some forms of cancer. The Body Mass Index (BMI) is the most widely accepted measurement of weight according to health professionals. Use the chart below to determine your own BMI by matching your height in the left hand column with your weight in the center. Anyone, male or female, with a BMI of 25 to 29 is considered over weight and a body mass index of 30 is considered obese. Remember to eat right and exercise.

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