Precautions

1. Clean exposed parts with a soft, slightly damp cloth. To remove stains, use a mild soap. Never use detergents, excess water, treated cloths, harsh cleaning agents, or sprays. Do not attempt to lubricate or open scale casing as this will void your warranty.
2. Although your scale is designed to be maintenance-free, it contains sensitive electronic parts. Avoid rough treatment. Do not drop, kick, or jump on it. Treat it with care to ensure the best performance.
3. This scale is designed to read weights between 25 and 440 lbs / 12 and 200 kgs. To prevent damage, do not place more than 440 lbs / 200 kg on it.
4. Do not store scale where you store cleaning chemicals. The vapors from some household products may affect the electronic components of your scale. Do not store scale on its side.
5. This scale is an extremely sensitive weighing device. In order to prevent run down of the batteries, do not store anything on the scale.
6. Do not dispose of batteries in fire. Batteries may explode or leak. Do not mix Alkaline, carbon-zinc (standard) or Nickel-Cadmium (rechargeable) batteries. Remove the batteries if the scale will not be used for a long period of time.

Five (5) Year Limited Warranty

This scale is warranted against defects in materials of workmanship (excluding batteries) for five (5) years of the original purchaser from date of retail purchase. It does not cover damages or wear resulting from accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair. Do not return to retailer. Should this scale require service (or replacement at our option) while under warranty, please pack the item in the original packaging and return it prepaid, along with store receipt showing date of purchase and a note explaining reason for return to:

TAYLOR Precision Products
2220 Entrada Del Sol, Suite A
Las Cruces, New Mexico 88001 USA

There are no express warranties except as listed above. This warranty gives you specific legal rights, and you may have other rights which vary from state to state. For additional product information, or warranty information outside the USA, please contact us through www.taylorusa.com.

©2011 Taylor Precision Products and its affiliated companies, all rights reserved. Taylor® and Leading the Way in Accuracy® are registered trademarks of Taylor Precision Products and its affiliated companies. All rights reserved.

Made to our exact specifications in China.

Not legal for trade.
www.taylorusa.com

CP7405 10.11

TAYLOR® 7405

Electronic Scale

Leading the Way in Accuracy®

Instruction Manual

Thank you for purchasing a Taylor® Electronic Scale. All of our precision instruments for weight measurement are engineered and designed to meet the highest quality standards...to assure you years of uncompromising accuracy and consistently dependable, convenient performance. Please read this instruction manual carefully before use. Keep these instructions handy for future reference.

General Operation

1. This scale operates on 3 "AAA" batteries, which are included. Please see "Battery Installation" section of this manual for installation instructions. Also, some models have a static cling label on the lens to prevent scratching. Please remove before use.
2. Your scale has been set at the factory to read weight in pounds. You may select kilogram of measurement by pressing the button on the bottom of the scale.
3. Place scale on flat, hard surface. Carpeted or uneven floors may affect accuracy.
4. Step on the scale (disregarding the reading) and step off. This process insures the scale is set and is now ready for use. Do this whenever the scale is moved. At all other times step straight onto the scale.
5. Step onto the scale platform. The display will count up to your weight. The display will flash twice when the weight is locked. Your weight is now displayed.
6. The scale will turn off automatically. **NOTE:** If you've weighed yourself and a "C" is displayed after the zero, weigh yourself again for a more accurate reading. If the scale displays "Err", this indicates the weight is not stable and the scale cannot lock in the weight. Step off and back onto the scale, standing still while your weight computes.

Low Battery

Your Taylor® scale is equipped with a low battery indicator. It will display "C" or readings will grow dim or become irregular if the batteries needs replacement.

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
Battery Installation

- Use 3 AAA batteries (included). Remove any plastic wrap from the batteries before proceeding. Always replace all 3 batteries at the same time; do not combine old and new batteries.
- Open the battery compartment on the bottom of the scale.
- Place the new batteries into the battery compartment according to the "+" and "-" polarity symbols at the bottom of the compartment.
- Replace the battery cover.


Trouble Shooting

In Case of Difficulty...

- Place your feet in the same position every time you weigh yourself. Stand still during measurement.
- Weigh yourself without clothing whenever possible. Clothes can vary in weight.
- Move the scale to a location where it won't be bumped. If your scale is bumped while being activated or while in use, an inaccurate reading may occur.
- For the most consistent readings, weigh yourself at the same time of day. Allow for weight fluctuations due to food or drink in your system.
- The maximum capacity of this scale is 440 lbs / 200 kgs. Should "O-LD" appear on the display during a weighing, this indicates that the maximum weighing capacity of the scale has been exceeded. Remove the weight immediately; otherwise, permanent damage may occur.

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
Reaching Target Weight

Maintaining proper weight is essential to healthy living. Obesity has been linked to a higher risk of heart disease, diabetes and some forms of cancer. The Body Mass Index (BMI) is the most widely accepted measurement of weight according to health professionals.

Use the chart below to determine your own BMI by matching your height in the left hand column with your weight in the center. Anyone, male or female, with a BMI of 25 to 29 is considered over weight and a body mass index of 30 is considered obese. Remember to eat right and exercise.

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