TENS Product Line-Up



TENS Unit Comparison Chart				
	CONSIGNATION OF CHARLES OF CHARLE	Long Long Long Long Long Long Long Long	CORROLL CORROL	
Model Number	Pocket Pain Pro™ PM400	Max Power Relief™ PM500	Total Power + Heat PM800	Avail ™ PM601
Pain Modes	3 Modes (Arm/Shoulder, Lower Back, Leg/Foot)	6 Modes (Arm, Lower Back, Leg/Foot, Joint, Shoulder)	9 Modes (Arm, Lower Back, Leg/ Foot, Joint, Shoulder)	6 Modes (Lower Back, Shoulder, Joint, Arm, Leg, Microcurrent)
Massage-like Modes	2 Modes (Knead, Steady)	3 Modes (Tap, Knead, Rub)	3 Modes (Tap, Knead, Rub)	4 Modes (Steady, Knead, Tap, Ac
Intensity Levels	10	15 (2x the power of the Pocket Pain Pro)	20	20
Channels	Single Channel	Single Channel	Single Channel	Dual Channel
Bluetooth App Controlled				$\sqrt{}$
Wireless				$\sqrt{}$
Heat			\checkmark	
Pads	1 Pair Long Life Pads™ (PMLLPAD)	1 Pair Long Life Pads™ (PMLLPAD)	2 Pair Gel Pads (PMGEL)	2 Wireless Pads (1PMWPAD - M, 1PMWPAD - L)
Screen		Easy-to-Read LCD	Easy-to-Read LCD	Smart phone
Maximum Therapy Duration	15 Min	15 Min	30 Min	60 Min
Warranty	Unit-1 year limited warranty Long Life Pads™-30 day limited warranty	Unit-1 year limited warranty Long Life Pads™-30 day limited warranty	Unit-1 year limited warranty Get Refills-30 day limited war- ranty	Unit-1 year limited warranty Wireless Pads-30 day limited warranty
Battery	2 AAA batteries included	2 AAA batteries included	Rechargeable Lithium Ion Battery	Rechargeable Lithium Ion Battery
Belt Clip		$\sqrt{}$		
QSG & Pad Placement Guide	$\sqrt{}$	$\sqrt{}$	\checkmark	Within App
Pad Holder	$\sqrt{}$	$\sqrt{}$	$\sqrt{}$	$\sqrt{}$
Travel Case				$\sqrt{}$