

# RECOVERY IS KEY

## **Polar Roll-On Massager** - GF-POL-PDQ - Display GF-POLROM - Single

- **Stainless Roller Ball stays ice cold**
- **Comfortable, soft rubber grip fits in your palm**
- **Smooth roll-on action**
- **Ice-pack-type, cold rolling massage**
- **Durable TPE and stainless components**
- **Highly portable and compact**
- **Product Size: 3" x 3"**
- **Downloadable Training Poster**



Smooth action,  
stainless steel  
polar ball swivels  
in all directions

Gray retaining ring  
unscrews to remove  
the ball for cleaning  
or freezing



PDQ display holds 9 units  
9.75" x 10.25" x 7.75" high

Place in freezer for  
approximately  
2 hours between uses

